

MC2F57 - July 1, 2016

Item # MC2F57 was discontinued on July 1, 2016. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

OPTICAL CHOPPER SYSTEM AND CHOPPER WHEELS

- ▶ **Crystal-Stabilized, Phase-Locked Feedback Loop Suppresses Low Frequency Drift and Pulse Jitter**
- ▶ **Manual, USB, and External Trigger Controller**
- ▶ **7 Blades Covering Frequencies from 1 Hz to 10 kHz**



MC1F2



MC1F15



MC1F30



MC1F60



MC1F100



MC2F57



MC1F10



MC2000

MC2000 Optical Chopper System
 Includes Controller,
 Chopper Head, 10-Slot Blade
 and Interface Cable

[Hide Overview](#)

OVERVIEW

Features

- Crystal-Stabilized, Phase-Locked Feedback Loop Suppresses Low Frequency Drift and Phase Jitter
- Harmonic, Subharmonic, and Fractional Harmonic Chopping with Sum and Difference Reference Outputs
- Microprocessor Controlled
- 2-Frequency Blades Available for Pump-Probe and Other Nonlinear Experiments
- Save and Recall User Setups in Non-Volatile RAM
- USB Interface
- Control Software Package Available



Screenshot of the English GUI

The MC2000 Optical Chopper, utilizing advanced features to meet the most demanding applications, is a rotary chopper that is composed of a control console, motor head assembly and slotted chopper wheel. The MC2000 uses a phase-locked loop (PLL) motor speed control design to precisely maintain the chopping speed and phase relative to a reference signal. An internal, crystal-stabilized, frequency synthesizer provides an accurate and stable reference frequency for stable long-term performance.

An optical switch is used to get the reference signal when the chopper blades pass through and interrupt the beam path. The switch incorporates

MC2000 Optical Chopper System Performance

| | |
|--|---|
| Chopping Frequency, with Various Blades | 1 Hz – 10 kHz |
| Frequency Drift | <20 ppm/°C |
| Ext. Reference Compatibility | TTL/CMOS |
| Frequency Resolution | 0.01 Hz (2-Slot Blade) 1 Hz (All Other Blades) |
| External Reference Signal Synchronization | |
| Harmonic | 2 to 15X |
| Subharmonic | 1/2 to 1/15X |

Chopper Wheel Options

| Item # | # of Slots | Chopping Frequency | Phase Jitter Max/Typ. ^c |
|---------------------|------------|--------------------|------------------------------------|
| MC1F2 ^a | 2 | 1 Hz – 99 Hz | 0.05° rms (Max) |
| MC1F10 ^b | 10 | 20 Hz – 1 kHz | 0.42°/0.13° rms |
| MC1F15 | 15 | 30 Hz – 1.5 kHz | 0.68°/0.27° rms |
| MC1F30 | 30 | 60 Hz – 3 kHz | 1.10°/0.45° rms |

an LED with a wavelength range of 850 nm - 940 nm that is operated at the 15 mA threshold current to minimize light output from the LED.

Unlike conventional, open-loop speed control designs, the PLL speed control circuit also allows the MC2000 chopper to be synchronized to external reference signals, including other MC2000 choppers and reference sources such as DSP lock-in amplifiers. Please see the *Operations* tab for more details on the MC2000 chopper and its applications.

The MC2000 design features a high quality, rare earth magnet, DC motor. The photo-etched optical chopper wheel blade is made of 0.01" thick blued clock spring steel for high precision. The compact optical head has a wide base for extra stability, and the base is slotted for two 1/4"-20 mounting screws on 2" centers. The interface cable uses circular snap-on LIMO connectors for easy setup.

The MC2000 controller includes a 240 x 128 pixel graphics display for setting and monitoring chopper functionality. All of the functions are accessible through a front panel control knob with turn and push control. Multiple user setups can be easily saved and recalled from non-volatile memory. A USB interface is included as a standard feature for remote PC control of the MC2000. A software package with a GUI to control the MC2000 can be downloaded from the *Software* tab. The package also supports controlling the MC2000 via programs written in C++ or LabVIEW.

Current owners of an MC2000 controller may upgrade the unit with the newest firmware release, which can be downloaded by clicking on the link on the *Software Downloads* tab. Please note that the system must remain powered on and in an idle state until the firmware update is complete.

Please Note: The MC2000 should be shielded from stray light sources because light contamination in the optical sensor will cause jitter.

| | | | |
|----------------|----------------------|--|-----------------|
| MC1F60 | 60 | 120 Hz – 6 kHz | 1.10°/0.78° rms |
| MC1F100 | 100 | 250 Hz – 10 kHz | 1.30°/1.10° rms |
| MC2F57 | Outer: 7 Inner: 5 | Outer: 14 – 700 Hz Inner: 10 – 500 Hz | 0.38°/0.08° rms |

- Please note the MC1F2 two-slot blade was updated to a newer version on September 21, 2011. This new blade is compatible with all units. If you purchased the older blade, please note that it can only be used with firmware version 3.11 and older. If you would like to update to the latest firmware and own the older style blade, please contact Tech Support for a free MC1F2 blade exchange.
- The MC1F10 blade is supplied with the unit. All other blades specified may be purchased separately through Thorlabs.
- Measured using internal reference signal.

[Hide Specs](#)

S P E C S

| Optical Head Specifications | |
|-----------------------------------|---|
| Chopping Blade Diameter | Ø4.0" (Ø101.6 mm) |
| Chopping Blade Thickness | 0.010" (0.254 mm) |
| Mounting Base | 1/4"-20 (or M8) Clearance Slots Spaced 3.0" (Compatible with Thorlabs Breadboards) |
| Mounting Hole | 1/4"-20 with 1/4" Max Screw Depth |
| Blade Specifications | |
| Chopping Blade Slots ^d | |
| MC1F2^a | 2 |
| MC1F10 (Default Blade) | 10 |
| MC1F15 | 15 |
| MC1F30 | 30 |
| MC1F60 | 60 |
| MC1F100 | 100 |
| MC2F57 | 7 Outer, 5 Inner |
| Slot Angle | |
| MC1F2 | 180° |
| MC1F10 (Default Blade) | 36° |
| MC1F15 | 24° |
| MC1F30 | 12° |
| MC1F60 | 6° |
| MC1F100 | 3.6° |
| MC2F57 | 51.4° Outer, 72° Inner |

| Physical Features | |
|-------------------------------------|--|
| Dimensions (W x H x D) | 5.8" x 2.8" x 12.5" (147 mm x 71 mm x 317.5 mm) |
| Input and Output Connectors | BNC |
| Menu Control | Twist / Push-Button Knob |
| Input Power Connection ^e | IEC Connector |
| Weight | 5 lbs (9.1 lbs Shipped Weight) |
| Operating Temperature | 10 – 40 °C |
| Display Type | 240 x 124 Pixel LCD Graphics Display |
| Frequency Resolution | 0.01 Hz (2 slot blade) 1 Hz (all other blades) |
| Power Supply | |
| Supply Type | Linear |
| Voltage Selection | 115 / 230 VAC Switchable |
| Input Voltage | 100/115 VAC ± 10%, 230 VAC ± 10% |
| Line Frequency | 50 – 60 Hz |
| Input Power | 20 VA Max |
| Fuse Ratings | 250 mA @ 115 VAC 125 mA @ 230 VAC |
| Fuse Type | IEC60127-2/III (250 V, Slo-blo Type 'T') |
| Fuse Size | 5 x 20 mm |

| Performance Specifications | |
|---------------------------------------|--|
| Chopping Frequency | |
| MC1F2 (2 slot) | 1 Hz – 99 Hz |
| MC1F10 (10 slot, Default Blade) | 20 Hz – 1 kHz |
| MC1F15 (15 slot) | 30 Hz – 1.5 kHz |
| MC1F30 (30 slot) | 60 Hz – 3 kHz |
| MC1F60 (60 slot) | 120 Hz – 6 kHz |
| MC1F100 (100 slot) | 250 Hz – 10 kHz |
| MC2F57 (2f slot) | Outer: 14 – 700 Hz Inner: 10 – 500 Hz |
| Phase Jitter (@ Max Frequency) | |
| MC1F2 (2 slot) ^a | 0.05° rms Max |
| MC1F10 (10 slot, Default Blade) | 0.42° rms (0.13° rms Typ.) |
| MC1F15 (15 slot) | 0.68° rms (0.27° rms Typ.) |
| MC1F30 (30 slot) | 1.10° rms (0.45° rms Typ.) |
| MC1F60 (60 slot) | 1.10° rms (0.78° rms Typ.) |
| MC1F100 (100 slot) | 1.30° rms (1.10° rms Typ.) |
| MC2F57 (2f slot) | 0.38° rms (0.08° rms Typ.) |
| Frequency Drift | <20 ppm/°C |
| Chopping Range | |
| Harmonic | 2 to 15x |
| Sub-Harmonic | 1/2 to 1/15x |
| Input/Output Specifications | |
| Ext. Input Compatibility | TTL/CMOS |
| Ext. Input Voltage Range ^b | 0 – 5 V |

| | |
|------------------------------------|--|
| Input High | >2 V |
| Input Low | <0.8 V |
| Ext. Input Impedance | 200 Ω |
| Ref Out Compatibility | TTL/CMOS |
| Ref Out Voltage Range ^b | 0 – 5 V Typical |
| Ref Out Impedance | 200 Ω |
| Min Load Impedance ^c | 500 Ω |
| Ref Out Signals | Inner/Outer Slot Chopping Blade, Synthesizer, Sum & Diff Frequencies |
| Ref Out Selection | Selectable Menu or USB command 'O' |
| Communications | |
| Communications Port | USB |
| Protocol | USB (RS232 Emulated) |
| Baud Rate | 115,200 (fixed) |
| Data Bits | 8 |
| Stop Bits | 1 |
| Parity | None |
| Handshaking | None |

- Please note the MC1F2 two-slot blade was updated to a newer version on September 21, 2011. This new blade is compatible with all units. If you purchased the older blade, please note that it can only be used with firmware version 3.11 and older. If you would like to update to the latest firmware and own the older style blade, please contact Tech Support for a free MC1F2 blade exchange.
- The reference output and external input is short circuit protected by limiting the current to 25 mA. Overvoltage and undervoltage protection is available, but continued use will degrade or damage the unit.
- The Min Load Impedance represents the smallest allowable terminating resistance. Applying lower impedances will cause the short circuit protection to limit the output voltage. Continued use in this mode will cause circuit degradation and eventual circuit failure.
- The MC1F10 blade is supplied with the unit. All other blades specified may be purchased separately through Thorlabs.
- The MC2000 is shipped with a 120 VAC US style power cord and the MC2000-EC is shipped with a 230 VAC power cord for use in Europe.

[Hide Pin Diagrams](#)

PIN DIAGRAMS

Controller Front Panel



[Click to Enlarge](#)

The front panel of the MC2000 offers two female BNC ports, one for an external input and one for the reference output. Pin diagrams for both ports are shown below.

Controller Back Panel



[Click to Enlarge](#)

The back panel of the MC2000 offers a USB 2.0 port and a male HR10-7R-6S interface connector. Pin diagrams for both connectors are shown below.

Chopper Back Panel




Click to Enlarge

The chopper head back panel offers a male HR10-7R-6S interface connector. The included HR10-7R-6S cable connects the chopper head to the back panel of the controller. The pin diagram for the connector is shown below.

External Input^a


BNC Female



a. TTL/CMOS, 200 Ω Input Impedance, 0 - 5 V Input Range

Reference Output^b

BNC Female



b. TTL/CMOS, 200 Ω Output Impedance, 0 - 5 V Input Range, 500 Ω Min Load Impedance

Computer Connection


USB Type B^c



c. USB 2.0 Type A to B Cable

Interface Connector

HR10-7R-6S




| Pin | Description | Pin | Description |
|-----|-------------|-----|-----------------|
| 1 | Motor + | 4 | Feedback Sensor |
| 2 | Motor - | 5 | Phase Sensor |
| 3 | VCC | 6 | Ground |

[Hide Operation](#)

OPERATION

Operation



An optical chopper is an electromechanical instrument that periodically interrupts a beam of light. The MC2000 Optical Chopper is a rotary type composed of a control console, motor head assembly, and slotted chopper wheel.

Applications

Unlike conventional, open-loop speed control designs, the PLL speed control circuit also allows the MC2000 chopper to be synchronized to external reference signals, including other MC2000 choppers and reference sources such as DSP lock-in amplifiers.

For more advanced measurements, the MC2000 can lock to a harmonic, sub-harmonic, or fractional-harmonic of an external reference frequency. The on board microprocessor is used to multiply the external reference up to the 15th harmonic or to divide the reference down to the 15th sub-harmonic. By combining both the frequency multiplication and division together, a fractional harmonic can be obtained.



MC2F57 2
Frequency Blade

In order to maximize the slope of the signal as the blade rotates (i.e., minimize the time to fully open or closed), the beam should be focused onto the widest portion of the slot in the chopper blade.

The MC2000 also supports 2-frequency chopping from a single chopper blade (shown at right). A special blade is available with 7 outer slots and 5 inner slots. This slot combination allows a single beam to be split and individually modulated for ratiometric experiments. Another example application is a pump-probe experiment where the pump beam is modulated at the outer frequency while modulating a probe beam at the inner frequency. The MC2000 provides the sum and difference frequencies of the 2-frequency blade for accurate lock-in detection of the frequency-mixed response.

[Hide Language Option](#)

LANGUAGE OPTION



Screenshot of the Chinese GUI

New Chinese Language Option

As part of our new multilingual interface initiative, we have released the MC2000 with a Chinese language display option. The display can be changed by sending a command using the command line language through the USB interface. The language display setting is stored in the unit's non-volatile memory. Current MC2000 controllers can be upgraded with the new display language by downloading the newest firmware release found on the *Software Downloads* tab.

[Hide Software Downloads](#)

SOFTWARE DOWNLOADS

Software with GUI for Remote Operation and Firmware Updates for the MC2000 Chopper Controller

Software

Version 1.0

Software package to control the MC2000. The installer includes the necessary drivers. This software package is compatible with the latest firmware for the system and supports writing programs to control the MC2000 using C++ or LabVIEW.

Firmware

Version 3.15/4.15

The MC2000 hardware has been revised several times and each revision requires different firmware versions. Hardware versions 3.0 and 4.0 are eligible for the firmware update available below.

To determine your hardware version, please attempt to update the firmware. Users of version 1.0 or 2.0 devices will be prompted that the firmware cannot be updated. All version 1.0 and 2.0 devices are out of warranty, and are no longer compatible with the latest

firmware.



[Hide Optical Chopper System](#)

Optical Chopper System



- ▶ Harmonic, Subharmonic, and Fractional Harmonic Chopping
- ▶ Sum and Difference Reference Outputs
- ▶ 10-Slot Blade Included (MC1F10, See Below)

Compatibility with Other Vendors

The MC2000 is also compatible with other manufacturer's blades. Select Scitec Instrument blades, including the 300D3/30, 300D8/60, 300D4/7, and 300D53/60 will work with this chopper, but due to slight variations in the design, not all functions will be available. Please note that we do not guarantee the functioning of other manufacturer's blades with our chopper head and damage may occur.

Existing MC2000 controllers can be upgraded by downloading the newest firmware release. Please note that the system must remain powered on and in an idle state until the firmware update is complete.

The MC2000 is shipped with a 120 VAC power cord for use in the United States while the MC2000-EC is shipped with a 230 VAC power cord for use in Europe. If you require a different power cord, please contact Tech Support prior to ordering.

| Part Number | Description | Price | Availability |
|-------------|--|------------|--------------|
| MC2000-EC | Optical Chopper System with MC1F10 10-slot (36°) Chopper Blade, 230 VAC Power Cord | \$1,200.00 | Today |
| MC2000 | Optical Chopper System with MC1F10 10-slot (36°) Chopper Blade, 120 VAC Power Cord | \$1,200.00 | Today |

[Hide Optical Chopper Wheels](#)

Optical Chopper Wheels

| | | | | | | | | | |
|--|------------------|-------------------|------------------|------------------|------------------|-----------------|-------------------|-------------------|--|
| Click to Enlarge | | | | | | | | | |
| Item # | MC1F2 | MC1F10 | MC1F15 | MC1F30 | MC1F60 | MC1F100 | MC2F57 | | |
| Slots | 2 | 10 | 15 | 30 | 60 | 100 | 7 Outer | 5 Inner | |
| Frequency | 1 - 99 Hz | 20 Hz - 1 kHz | 30 Hz - 1.5 kHz | 60 Hz - 3 kHz | 120 Hz - 6 kHz | 250 Hz - 10 kHz | 2 Frequencies | | |
| Laser Beam Diameter^a | | | | | | | | | |
| Max. Possible ^b | | 13.23 mm (0.521") | 9.25 mm (0.364") | 4.86 mm (0.191") | 2.49 mm (0.098") | 1.51 mm (0.06") | | | |
| Max. Center ^c | 26.1 mm (1.028") | 10.32 mm (0.406") | 6.9 mm (0.272") | 3.45 mm (0.136") | 1.73 mm (0.068") | 1.04 mm (0.04") | 14.96 mm (0.589") | 14.35 mm (0.565") | |
| Max. Slot ^d | | 6.34 mm (0.249") | 3.99 mm (0.157") | 1.89 mm (0.074") | 0.92 mm (0.036") | 0.55 mm (0.02") | | | |

- In order to maximize the slope of the signal as the blade rotates (i.e., minimize the time to fully open or closed), the beam should be focused onto the widest portion of the slot.
- The largest beam diameter that can pass through the widest portion of the chopper wheel slot, indicated by the red circle in the drawing to the right.
- The largest beam diameter that can pass through the center of the chopper wheel slot, indicated by the green circle in the drawing to the right.



Click to Enlarge

- The largest beam diameter that can pass through the narrowest portion of the chopper wheel slot, indicated by the blue circle in the drawing to the right.

Thorlabs offers a number of different blades to extend the range of the chopping frequency. A standard 10-slot blade (MC1F10) is included with the MC2000. 2-slot, 15-slot, 30-slot, 60-slot, 100-slot, and 2-frequency blades are also available. Our blades are manufactured from 0.01" thick blued clock spring steel.

The MC2000 also supports 2-frequency chopping from a single chopper blade. A special blade is available with 7 outer slots and 5 inner slots. This slot combination allows a single beam to be split and individually modulated for ratiometric experiments. Please see the *Operation* tab for more details.

Please note the MC1F2 two-slot blade was updated to a newer version on September 21, 2011. This new blade is compatible with all units. If you purchased the older blade, please note that it can only be used with firmware version 3.11 and older. If you would like to update to the latest firmware and own the older style blade, please contact Tech Support for a free MC1F2 blade exchange.

| Part Number | Description | Price | Availability |
|-------------|---|---------|--------------|
| MC1F2 | 2 Slot Blade for Optical Chopper, 1 to 99 Hz | \$58.90 | Today |
| MC1F10 | 10 Slot Blade for Optical Chopper, 20 to 1,000 Hz | \$41.80 | Today |
| MC1F15 | 15 Slot Blade for Optical Chopper, 30 to 1,500 Hz | \$41.80 | Today |
| MC1F30 | 30 Slot Blade for Optical Chopper, 60 to 3,000 Hz | \$41.80 | Today |
| MC1F60 | 60 Slot Blade for Optical Chopper, 120 to 6,000 Hz | \$41.80 | Today |
| MC1F100 | 100 Slot Blade for Optical Chopper, 250 to 10,000 Hz | \$45.80 | Today |
| MC2F57 | 2 Frequency, 7 Outer/5 Inner Slot Blade for Optical Chopper | \$60.00 | Today |

Visit the *Optical Chopper System and Chopper Wheels* page for pricing and availability information:
https://www.thorlabs.com/newgrouppage9.cfm?objectgroup_id=287