


## Legacy LEMO Stepper Motor Drive Cables

The conversion cable, PAA621, can be used to connect a 15-pin D-type connector to the 8-pin LEMO connector on any of Thorlabs' stepper motor actuators.



PAA621



ITEM#	\$	£	€	RMB	DESCRIPTION
PAA621	\$ 96.95	£ 67.30	€ 86,10	¥ 818.70	Stepper Motor Drive Cable, Female D-Type to Male LEMO Converter, 3.0 m

## Legacy 25-Pin D-Type Stepper Motor Cable

The PAA611 cable features a 25-pin male D-type connector on one end and a 15-pin female D-type connector on the other end. This cable can be used with the MST601, BSC001, BSC002, and BSC003 stepper motor controllers.



PAA611



ITEM#	\$	£	€	RMB	DESCRIPTION
PAA611	\$ 61.80	£ 42.90	€ 54,90	¥ 521.90	Legacy APT Stepper Motor Drive Cable, 3.0 m

## Z8 DC Motor Extension Cable

The PAA632 is an extension cable for use with the Z8 series actuators (page 534). This cable is 2.5 m (8.2') long with a 15-pin D-type male connector on one end and a female 15-pin D-type connector on the other end.



PAA632

NEW  
product

ITEM#	\$	£	€	RMB	DESCRIPTION
PAA632	\$ 52.80	£ 36.70	€ 46,90	¥ 445.90	Z8 Series Extension Cable, 2.5 m

## Laser Diode and Temperature Controller Cables

These cables are the same cables included with many of our laser diode current and TEC controllers. The table below outlines the controller that each cable is typically used with.

ITEM#	USE	COMPATIBLE CONTROLLER
CAB400	LD	LDC200C Series, PRO 8 Series
CAB420-15	TEC	TED200C, TTC001, PRO8 Series
CAB430	LD and TEC	ITC100 Series
CAB530	LD and TEC	ITC5022, ITC5052, ITC5102
CAB4000	TEC	TED4015, ITC4000 Series – 5 A Max
CAB4001	TEC	TED4015 ITC4000 Series – 20 A Max
CAB4005	LD	LDC4000 Series, ITC4000 Series – 5 A Max
CAB4006	LD	LDC4000 Series, ITC4000 Series – 20 A Max

CAB420-15

NEW  
versions

CAB4001

ITEM#	\$	£	€	RMB	DESCRIPTION
CAB400	\$ 66.00	£ 45.80	€ 58,60	¥ 557.40	LD Current Controller Cable, 9-Pin D-Sub, 1.5 m
CAB420-15	\$ 72.00	£ 50.00	€ 64,00	¥ 608.00	TEC Controller Cable, 15-Pin D-Sub, 1.5 m
CAB430	\$ 120.00	£ 83.20	€ 106,60	¥ 1,013.30	LD and TEC Controller Cable for ITC100, 1.5 m
CAB530	\$ 96.00	£ 66.60	€ 85,30	¥ 810.70	ITC5000 Series to TEC & LD Driver Mounts
CAB4000	\$ 65.00	£ 45.10	€ 57,80	¥ 548.90	TEC Controller Cable for TED4000, 5 A Max, 9-Pin D-Sub, 1.5 m
CAB4001	\$ 170.00	£ 117.90	€ 151,00	¥ 1,435.50	TEC Controller Cable for TED4000, 20 A Max, 17W2, 1.5 m
CAB4005	\$ 80.65	£ 56.00	€ 71,70	¥ 681.10	Cable for LDC4000 Series, 5 A, 13W3 to D-Sub-9, 1.5 m
CAB4006	\$ 80.65	£ 56.00	€ 71,70	¥ 681.10	Cable for LDC4000 Series, 20 A, 13W3 to 13W3, 1.5 m