

**Performance™ Series Optical Breadboards – 60mm Thick**

**Metric Performance™ Breadboards, 60mm Thick – Stocked for Immediate Shipment**

| ITEM#    | \$          | £          | €          | RMB         | DIMENSIONS         | WEIGHT       |
|----------|-------------|------------|------------|-------------|--------------------|--------------|
| PBH51501 | \$ 409.70   | £ 258.10   | € 381.00   | ¥ 3,912.60  | 300 x 300 x 60mm   | 24lbs/11kg   |
| PBH51502 | \$ 474.30   | £ 298.80   | € 441.10   | ¥ 4,529.60  | 600 x 300 x 60mm   | 44lbs/20kg   |
| PBH51528 | \$ 559.30   | £ 352.40   | € 520.10   | ¥ 5,341.30  | 600 x 450 x 60mm   | 75lbs/34kg   |
| PBH51505 | \$ 635.80   | £ 400.60   | € 591.30   | ¥ 6,071.90  | 600 x 600 x 60mm   | 82lbs/37kg   |
| PBH51510 | \$ 793.90   | £ 500.20   | € 738.30   | ¥ 7,581.70  | 750 x 750 x 60mm   | 123lbs/56kg  |
| PBH51506 | \$ 775.20   | £ 488.40   | € 720.90   | ¥ 7,403.20  | 900 x 600 x 60mm   | 119lbs/54kg  |
| PBH51511 | \$ 1,028.50 | £ 648.00   | € 956.50   | ¥ 9,822.20  | 900 x 750 x 60mm   | 148lbs/67kg  |
| PBH51514 | \$ 1,040.40 | £ 655.50   | € 967.60   | ¥ 9,935.80  | 900 x 900 x 60mm   | 174lbs/79kg  |
| PBH51507 | \$ 936.70   | £ 590.10   | € 871.10   | ¥ 8,945.50  | 1200 x 600 x 60mm  | 157lbs/71kg  |
| PBH51512 | \$ 1,139.00 | £ 717.60   | € 1,059.30 | ¥ 10,877.50 | 1200 x 750 x 60mm  | 194lbs/88kg  |
| PBH51515 | \$ 1,229.10 | £ 774.30   | € 1,143.10 | ¥ 11,737.90 | 1200 x 900 x 60mm  | 231lbs/105kg |
| PBH51508 | \$ 1,059.10 | £ 667.20   | € 985.00   | ¥ 10,114.40 | 1500 x 600 x 60mm  | 194lbs/88kg  |
| PBH51513 | \$ 1,336.20 | £ 841.80   | € 1,242.70 | ¥ 12,760.70 | 1500 x 750 x 60mm  | 240lbs/109kg |
| PBH51516 | \$ 1,429.70 | £ 900.70   | € 1,329.60 | ¥ 13,653.60 | 1500 x 900 x 60mm  | 265lbs/120kg |
| PBH51517 | \$ 1,676.20 | £ 1,056.00 | € 1,558.90 | ¥ 16,007.70 | 1800 x 900 x 60mm  | 342lbs/155kg |
| PBH51525 | \$ 2,108.00 | £ 1,328.00 | € 1,960.40 | ¥ 20,131.40 | 1800 x 1250 x 60mm | 472lbs/214kg |
| PBH51522 | \$ 1,929.50 | £ 1,215.60 | € 1,794.40 | ¥ 18,426.70 | 2000 x 1000 x 60mm | 485lbs/220kg |
| PBH51503 | \$ 584.80   | £ 368.40   | € 543.90   | ¥ 5,584.80  | 900 x 300 x 60mm   | 62lbs/28kg   |
| PBH51504 | \$ 673.20   | £ 424.10   | € 626.10   | ¥ 6,429.10  | 1200 x 300 x 60mm  | 82lbs/37kg   |
| PBH51521 | \$ 1,608.20 | £ 1,013.20 | € 1,495.60 | ¥ 15,358.30 | 1500 x 1000 x 60mm | 317lbs/144kg |
| PBH51509 | \$ 1,186.60 | £ 747.60   | € 1,103.50 | ¥ 11,332.00 | 1800 x 600 x 60mm  | 231lbs/105kg |

\*These breadboards can be used with our ScienceDesk™ workstations. See pages 2-7 to learn more.

ScienceDesk™

Breadboards

Breadboard Supports

Optical Tables

Table Supports

Non-Magnetic

**Expanding The Line**



Our ScienceDesk™ frames are designed to accept breadboards from our UltraLight, Performance, and PerformancePlus lines, thus allowing flexibility in creating a workstation with an array of vibration isolation possibilities.

- Wide Range of Accessories
- Six Standard Frame Sizes With 765mm (30.6") Height
  - 600 x 900mm
  - 750 x 750mm
  - 750 x 900mm
  - 750 x 1200mm
  - 900 x 900mm
  - 900 x 1200mm
- Active, Passive, and Rigid Isolation Systems
- Made From ESD Materials
- Ergonomic Design
- Heavy-Duty Casters

**ScienceDesk™ and Accessories**



PC Moveable Shelf



Side Shelf



Under Shelf



Single Shelf



Adjustable Wrist Rest

There is an array of accessories for the ScienceDesk™ workstations including under shelves, moveable shelves, and ergonomic accessories such as wrist rests and keyboard holders.

**See Pages 1-14 for More Information**