

Mounting Adapter for MAX200 Series Stages for Nikon TE2000 Microscope


Features

- Pre-Drilled Tap Holes for Easy Mounting of Motorized Sample Stage
- Mounts MAX200 Stages to the Nikon TE2000 Microscope




The MAX200A1 is an adapter for easy mounting of our MAX200 series stages to the Nikon TE2000 microscope.

ITEM#	\$	£	€	RMB	DESCRIPTION
MAX200A1	\$ 125.00	£ 78.80	€ 116,30	¥ 1,193.80	Mounting Adapter for Nikon TE2000 Microscope




MAX200P1 Features

- Compatible With the MAX200 and MAX201 Motorized Stages
- Fits Standard Microscope Slides and Ø35mm Petri Dishes
- Clip Holder to Secure Samples
- Size: 128mm x 136mm x 2mm (5.04" x 5.35" x 0.08")



MAX200P6 Features

- Compatible With the MAX202 and MAX203 Motorized Stages
- Fits Standard Microscope Slides and Ø35mm Petri Dishes
- Clip Holder to Secure Samples
- Size: 86mm x 143mm x 2mm (3.4" x 5.65" x 0.08")



MAX200P2 Features

- Compatible With the MAX200 and MAX201 Motorized Stages
- Accepts Standard Microscope Slides and Ø35mm Petri Dishes
- Tapped Holes: 63 x 1/4"-20 (M6) and 86 #4-40 (M3) Tapped Holes, 1/2" (12.5mm) Pitch
- Size: 129mm x 137mm x 6mm (5.1" x 5.4" x 0.24")



MAX200P3 Features

- Compatible With the MAX200 and MAX201 Motorized Stages
- Accepts Standard Microscope Slides
- Compatible With Ø35mm Petri Dishes
- Magnetic Steel Platform
- Size: 129 x 137 x 2mm (5.1" x 5.4" x 0.08")

ITEM#	\$	£	€	RMB	DESCRIPTION	COMPATIBILITY
MAX200P1	\$ 200.00	£ 126.00	€ 186.00	¥ 1,910.00	Microscope Slide Platform for MAX200 & MAX201 Stages	MAX200 & MAX201
MAX200P6	\$ 200.00	£ 126.00	€ 186.00	¥ 1,910.00	Microscope Slide Platform for MAX202 & MAX203 Stages	MAX202 & MAX203
MAX200P2	\$ 125.00	£ 78.80	€ 116.30	¥ 1,193.80	Tapped Hole Microscope Sample Platform (1/4-20 & 4-40)	MAX200 & MAX201
MAX200P2/M	\$ 125.00	£ 78.80	€ 116.30	¥ 1,193.80	Tapped Hole Microscope Sample Platform (M6 & M3)	MAX200 & MAX201
MAX200P3	\$ 95.00	£ 59.90	€ 88.40	¥ 907.30	Magnetic Microscope Sample Platform	MAX200 & MAX201

Positioning Component Compatibility

POSITIONING STAGE	OPTIONAL PIEZO STAGES	COMPATIBLE SAMPLE TYPE/PLATFORM	Olympus IX71	Zeiss Axiovert 200	Nikon TE2000
MAX202	N/A	Well Plate: No Platform Required Slide or Petri Dish Holder: MAX200P6	yes	yes	yes*
MAX203	N/A	Well Plate: No Platform Required Slide or Petri Dish Holder: MAX200P6	yes	yes	yes*
MAX200 MAX201	MicroScan Series Stages SCXY100, SCZ100, SCXYZ100 (See Page 630)	Slide or Petri Dish Holders: MAX200P1 Sample Holder MAX200P2 Threaded Sample Platform MAX200P3 Magnetic Sample Platform	yes	no	yes*

* The MAX200A1 is needed to mount the MAX200 series XY stepper stages to the Nikon TE2000 microscope.