

MPM-XYRS - April 18, 2019

Item # MPM-XYRS was discontinued on April 18, 2019. For informational purposes, this is a copy of the website content at that time and is valid only for the stated product.

MANUAL XY STAGE FOR THE NIKON FN1 MICROSCOPE

- ▶ Compatible with All Nikon FN1 Upright Microscopes
- ▶ Manual Adjustment in X and Y

Application Idea
MPM-XYRS Mounted
on a Nikon FN1 Microscope



MPM-XYRS



[Hide Overview](#)

OVERVIEW

The MPM-XYRS stage is designed as a manual XY stage for the Nikon FN1 microscope and is meant for use in physiology applications. It features an M6-tapped hole matrix on 25 mm centers on either side of a central sample plate. The stage features tension adjusters with 30 mm stroke in X and Y.

Thorlabs also offers a line of motorized microscopy stages for both XY and Z sample positioning.

[Hide Part Numbers](#)

Part Number	Description	Price	Availability
MPM-XYRS	Manual XY Stage for Nikon FN1 Microscopes	\$3,057.25	Lead Time

