

For current pricing, please see our website.

CHAPTERS

- Laser Scanning Microscopy
- Microscopy Components
- OCT Imaging Systems
- OCT Components
- Adaptive Optics


SECTIONS

- Microscopy Stages
- ScienceDesk
- LEDs
- Light Sources
- Objectives/Scan Lenses
- Dispersion Compensating Mirrors
- Fluorescence Imaging Filters
- Filter Cubes
- Scanning Mirrors
- PMT Modules
- Microscope Adapters
- Cuvette Holder
- FiberPorts
- Test Targets/Reticles
- ThorScopes
- Focus Blocks
- Pinhole Wheel

Mounting Brackets/Adapters for XY Microscopy Stages

These mounting brackets and adapters allow the MLS203-1 XY Microscopy Stage to be mounted to various upright and inverted microscopes from Nikon and Olympus. Please refer to the table below for specific information on compatibility. The MLS203-2 XY microscopy stage is directly compatible with Zeiss Axio Observer and Axiovert 40 microscopes and does not require the purchase of additional mounting brackets.

NEW products




The mounting brackets attach to the MLS203 stage using M5 x 12 bolts. Mounting of the MLSA02 is shown above.

MICROSCOPE TYPE	UPRIGHT/INVERTED	MOUNTING BRACKET
Olympus IX71 and IX81	Inverted	MLSA02
Olympus BX51 and BX61	Upright	MLSA05
Nikon TE2000 and Eclipse Ti	Inverted	MLSA03
Nikon 80i and 90i	Upright	MLSA06

Mounting Brackets for the MLS203-1 Series Stage

ITEM #	\$	£	€	RMB	DESCRIPTION
MLSA02	\$ 105.00	£ 75.60	€ 91,35	¥ 836.85	Olympus IX71/IX81 Mounting Bracket for MLS203-1 Stage
MLSA03	\$ 129.00	£ 92.88	€ 112,23	¥ 1,028.13	Nikon TE2000/Eclipse Ti Mounting Bracket for MLS203-1 Stage
MLSA05	\$ 350.00	£ 252.00	€ 304,50	¥ 2,789.50	Olympus BX51/61 Adapter for MLS203-1 Stage
MLSA06	\$ 350.00	£ 252.00	€ 304,50	¥ 2,789.50	Nikon 50/80/90i Adapter for MLS203-1 Stage

Breadboard Mounting Brackets



MLS203 Stage shown with MLSA01 Work Surface Adapters

The MLSA01 Brackets enable the MLS203-1 XY Microscopy Stage to be bolted to an optical table or breadboard, thereby facilitating the use of this stage in custom-built microscope setups or general photonics applications. The brackets, which are sold as a set of two, fit both imperial and metric mounting surfaces. Please note that these brackets are not compatible with the MLS203-2 Microscopy Stage.

NEW products



MLSA01
Breadboard Mounting Brackets

Breadboard Mounting Brackets

ITEM #	\$	£	€	RMB	DESCRIPTION
MLSA01	\$ 155.00	£ 111.60	€ 134,85	¥ 1,235.35	Breadboard/Optical Table Mounting Brackets