

SM2 Series Endcaps and End Plug

These endcaps and end plug have been designed to protect delicate optical assemblies from dust when not in use. The blank plug (SM2PL) contains a centered drill point to facilitate custom machining.



SM2CP1



SM2CP2



SM2PL

NEW
product




ITEM#	\$	£	€	RMB	DESCRIPTION
SM2CP1	\$ 17.75	£ 12.40	€ 15,80	¥ 149.90	SM2 Series Endcap, Internal Threads
SM2CP2	\$ 17.75	£ 12.40	€ 15,80	¥ 149.90	SM2 Series Endcap, External Threads
SM2PL	\$ 14.00	£ 9.70	€ 12,50	¥ 118.30	SM2 Series Machinable End Plug, External Threads

SM2 Series Lens Tube Couplers

- SM2TI: Connect Two Externally Threaded SM2 Components
- SM2T2: Connect Two Internally Threaded SM2 Components
- SM2NT: Limit the External Thread Depth of an SM2 Component

Also Available
in Ø1/2",
Ø1.00", and Ø3.00"

Please refer to our website for complete models and drawings.




ITEM#	\$	£	€	RMB	D	DESCRIPTION
SM2T1	\$ 16.75	£ 11.70	€ 14,90	¥ 141.50	Ø02.20"	SM2 Coupler, Internal Threads
SM2T2	\$ 33.70	£ 23.40	€ 30,00	¥ 284.60	Ø02.25"	SM2 Coupler, External Threads
SM2NT	\$ 13.30	£ 9.20	€ 11,90	¥ 112.40	Ø02.25"	SM2 (Ø2.035"-40) Threaded Locking Nut


SM2-Compatible Quick-Release Lens Tube Couplers



SM2Q

Easily Swap
ComponentsSM2QF
HousingSM2QM
Ring

Also Available in
Ø1/2", Ø1.00",
Ø30 mm, Ø3.00"

Quarter-Turn
Connection and Removal

- Quarter-Turn Connect and Disconnect Adapters
- Lockable Design
- Easily Swap SM2-Threaded Components

These adapters provide a means to quickly connect and disconnect two similarly SM2-threaded components such as lens tubes, mounted optics, and detectors. Thread a male adapter onto an externally SM2-threaded component and thread the female adapter onto an internally SM2-threaded component. A witness line on each adapter allows them to be aligned easily.

Please refer to our website for complete models and drawings.

ITEM#	\$	£	€	RMB	D	DESCRIPTION
SM2Q	\$ 43.20	£ 30.00	€ 38,40	¥ 364.80	Ø2.50" (Ø63.5 mm)	SM2 Quick-Release Adapter Set
SM2QM	\$ 24.00	£ 16.70	€ 21,40	¥ 202.70	—	SM2 Quick-Release Adapter, Male
SM2QF	\$ 24.00	£ 16.70	€ 21,40	¥ 202.70	Ø2.50" (Ø63.5 mm)	SM2 Quick-Release Adapter, Female