


**Cage  
Products  
Available  
in 16, 30,  
& 60 mm  
Versions**

## 30 mm and 60 mm Cage System Mounting Platform




- Large Mounting Surface with a Multitude of Various Mounting Holes
- Accepts Our Small and Large Clamping Arms (PM1 and PM2)
- Fits Both 30 mm and 60 mm Cage Systems




ITEM#	METRIC ITEM#	\$	£	€	RMB	DESCRIPTION
LCPM	LCPM/M	\$ 85.10	£ 59.00	€ 75,60	£ 718.60	Horizontal Cage System Mounting Platform
PM1	PM1/M	\$ 9.50	£ 6.60	€ 8,50	£ 80.30	Small Clamping Arm
PM2	PM2/M	\$ 14.70	£ 10.20	€ 13,10	£ 124.20	Large Clamping Arm

## 30 mm Cage System U-Bench

The CB1 Cage System U-Bench provides a rigid construction gap within a standard 30 mm Cage System. A generous supply of mounting holes are provided so that small optical assemblies can be constructed within the U-Bench that might not otherwise fit within the standard cage system.



- SM1 (Ø1.035"-40) Tapped Holes at Both Interface Ends
- Compatible with Standard 30 mm Cage Systems
- Solid, One-Piece Construction Maintains System Alignment



ITEM#	METRIC ITEM#	\$	£	€	RMB	DESCRIPTION
CB1	CB1/M	\$ 72.50	£ 50.30	€ 64,40	¥ 612.20	30 mm Cage System U-Bench

## 30 mm Cage-Compatible Photomultiplier Tube Adapter

**NEW  
product**



- Light-Tight Design
- TR and P-Post Mountable
- Securely Holds a PMT via Threaded Retaining Collar
- 30 mm Cage System Compatible

The PXT1 has been specifically designed to provide a light-tight mounting solution for integrating Ø1-1/8" barrel side-on photomultiplier tubes (PMTs) into Thorlabs' SM1 (Ø1.035"-40) Lens Tubes and 30 mm Cage Systems. A complete list of compatible Ø1-1/8" barrel, side-on Hamamatsu PMTs has been provided.

### Compatible Hamamatsu PMTs

- |          |          |           |           |
|----------|----------|-----------|-----------|
| ■ 1P21   | ■ R2658  | ■ R508    | ■ R7518   |
| ■ 1P21P  | ■ R2658R | ■ R5108   | ■ R7518P  |
| ■ 1P28   | ■ R2693  | ■ R5983   | ■ R758-10 |
| ■ 1P28A  | ■ R2693P | ■ R5983P  | ■ R7639   |
| ■ 931A   | ■ R2949  | ■ R5984   | ■ R8486   |
| ■ 931B   | ■ R3788  | ■ R636-10 | ■ R8487   |
| ■ R105   | ■ R3896  | ■ R636P   | ■ R9110   |
| ■ R105UH | ■ R4220  | ■ R6925   | ■ R9220   |
| ■ R1527  | ■ R4220P | ■ R7154   | ■ R928    |
| ■ R1527P | ■ R4332  | ■ R7447   | ■ R928P   |
| ■ R212   | ■ R446   | ■ R7447P  | ■ R955    |
| ■ R212UH | ■ R4632  | ■ R7517   | ■ R955P   |

ITEM#	\$	£	€	RMB	DESCRIPTION
PXT1	\$ 96.00	£ 66.60	€ 85,30	¥ 810.70	30 mm Cage-Compatible Photomultiplier Tube Adapter