


Free-Standing Instrumentation Shelf

The Thorlabs free-standing instrumentation shelf provides storage space or further mounting surfaces for instruments, power supplies, and other accessories. The system is available in two sizes that span 6.0' (2.0 m) and 8.0' (2.5 m). Simply match the shelf length to the length or width of your optical table. Each system comes complete with two end frames, two adjustable upper shelves, and one adjustable under table shelf, and a power distribution outlet.

The 12-outlet power strip is available in several power supply outlet types. Shelves fitted with a NEMA-type outlet (PTA261 and PTA271) are shipped from stock. Shelves fitted with IEC (PTA262 and PTA272), SCHUKO (PTA263 and PTA273), or UK (PTA264 and PTA274) type outlets are built upon receipt of order.

Features


- Stable, Free-Standing Design Does Not Interfere with Vibration Isolation
- Shelf Height is Adjustable from 43.21" to 70.79" (1.098 to 1.798 m)
- Power Strips Included (Choose from NEMA, IEC, SCHUKO, or UK)
- Both Overhead and Under Table Shelves Included
- **Frame Height:** 72.64" (1.845 m)
- **Shelf Height (Max):** 70.79" (1.798 m)
- **Load Capacity (Max):** 154 lbs (70 kg) per Shelf
- **Shelf Width:** 15.75" (400 mm)



PTA261



Power Strip Included
Four Options Available



Free-Standing Shelf

ITEM#	\$	£	€	RMB	DESCRIPTION
PTA261	\$ 1,385.00	£ 960.20	€ 1,229.70	¥ 11,695.00	Overhead Shelf, 2 x 0.4 m (NEMA)
PTA271	\$ 1,465.00	£ 1,015.50	€ 1,300.50	¥ 12,371.00	Overhead Shelf, 2.5 x 0.4 m (NEMA)

Stocked for Immediate Delivery

Free-Standing Shelf

ITEM#	\$	£	€	RMB	DESCRIPTION
PTA262	\$ 1,385.00	£ 960.20	€ 1,229.70	¥ 11,695.00	Overhead Shelf, 2 x 0.4 m (IEC)
PTA263	\$ 1,385.00	£ 960.20	€ 1,229.70	¥ 11,695.00	Overhead Shelf, 2 x 0.4 m (SCHUKO)
PTA264	\$ 1,385.00	£ 960.20	€ 1,229.70	¥ 11,695.00	Overhead Shelf, 2 x 0.4 m (UK)
PTA272	\$ 1,465.00	£ 1,015.50	€ 1,300.50	¥ 12,371.00	Overhead Shelf, 2.5 x 0.4 m (IEC)
PTA273	\$ 1,465.00	£ 1,015.50	€ 1,300.50	¥ 12,371.00	Overhead Shelf, 2.5 x 0.4 m (SCHUKO)
PTA274	\$ 1,465.00	£ 1,015.50	€ 1,300.50	¥ 12,371.00	Overhead Shelf, 2.5 x 0.4 m (UK)

Call for Lead Time