


Vega Series OCT Systems

VEGA SERIES PRECONFIGURED SYSTEM	
ITEM #	VEG110C1
Base Unit	VEG110
Center Wavelength	1300 nm
Scanning System	OCTG-1300NR (Standard Scanner)
Scan Lens Kit	OCT-LK4
Accessories	OCT-STAND(/M) and OCT-XYR1(/M)



VEG110C1
The VEG110C1 configuration consists of a VEG110 long imaging range base unit, standard scanner, long depth of focus scan lens kit, stand, and stage.

VEGA SERIES OCT BASE UNIT SPECIFICATIONS	
ITEM #	VEG110
Center Wavelength	1300 nm
Imaging Depth (Air/Water)	12 mm / 9.0 mm
Axial/Depth Resolution (Air/Water)	16 μ m / 12 μ m
A-Scan/Line Rate	100 kHz
Sensitivity*	101 dB (at 100 kHz)
Base Unit Dimensions	321 mm x 320 mm x 213 mm

*Typical Values Measured Using a Scanner With Dual Path Setup

