













A. TOP PLATE DETAIL





FOR INFORMATION ONLY
NOT FOR MANUFACTURING PURPOSES

DRAWING PROJECTION		
	NAME	DATE
DRAWN	CJL	04/AUG/11
APPROVAL	AR	04/AUG/11

COPYRIGHT © 2011 BY THORLABS

VALUES IN PARENTHESIS ARE CALCULATED AND MAY CONTAIN ROUNDOFF ERRORS

THORLABS
www.thorlabs.com

700mm (27.5") RIGID SUPPORT
FRAME 1500 x 750mm (5' x 2.5')

MATERIAL	REV
STEEL	A

ITEM #	APPROX WEIGHT
PFR75150-7	65.5 kg